

DENTAL
X-RAY EXPOSURE TECHNIQUES

<u>Exam</u>	<u>kVp</u>	<u>mA</u>	<u>Time</u>
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____

EMERGENCY AND SAFETY PROCEDURES

1. Should this machine malfunction immediately switch off the electrical power.
2. Move the patient if conditions make it advisable.
3. Notify the safety officer.

DENTAL
X-RAY EXPOSURE TECHNIQUES

<u>Exam</u>	<u>kVp</u>	<u>mA</u>	<u>Time</u>
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____

EMERGENCY AND SAFETY PROCEDURES

1. Should this machine malfunction immediately switch off the electrical power.
2. Move the patient if conditions make it advisable.
3. Notify the safety officer.